


Estimation of Ferrous Iron in Mohr's salt

Theory

Potassium dichromate quantitatively oxidises ferrous iron (Fe^{2+}) of Mohr's salt to ferric iron (Fe^{3+}) in acid medium.



The solution of Mohr's salt can be titrated with standard $\text{K}_2\text{Cr}_2\text{O}_7$ solution in the presence of H_3PO_4 and $2(\text{N})\text{H}_2\text{SO}_4$ using Barium diphenylamine sulfonate as indicator.

Procedure

1) Preparation of standard $(\frac{\text{N}}{10})\text{K}_2\text{Cr}_2\text{O}_7$ solution

Weigh out accurately 1.225 gm $\text{K}_2\text{Cr}_2\text{O}_7$ in a 250 ml volumetric flask. Dissolve it and make the volume up to the mark with distilled water.

Equivalent weight of $\text{K}_2\text{Cr}_2\text{O}_7$ solution = 49.03

1000 ml 1(N) $\text{K}_2\text{Cr}_2\text{O}_7$ solution = 49.03 gm of $\text{K}_2\text{Cr}_2\text{O}_7$

250 ml $(\frac{\text{N}}{10})\text{K}_2\text{Cr}_2\text{O}_7$ solution = 1.225 gm of $\text{K}_2\text{Cr}_2\text{O}_7$

Table 1: Preparation of standard ($N/10$) $K_2Cr_2O_7$ solution

Weight of the bottle with $K_2Cr_2O_7$ (W ₁ gm)	Weight of the bottle after transfer (W ₂ gm)	Weight taken (W ₁ - W ₂ = W gm)
.....

$$\text{Strength of } K_2Cr_2O_7 = \frac{\text{Weight taken}}{\text{Weight required}} \left(\frac{N}{10} \right)$$

$$= \frac{W_1 - W_2}{1.225} \left(\frac{N}{10} \right)$$

$$= \frac{W}{1.225} \left(\frac{N}{10} \right)$$

$$= f \left(\frac{N}{10} \right), \quad f = \text{factor}$$

2) Titration of Mohr's salt (Fe^{2+} solution) with standard $\text{K}_2\text{Cr}_2\text{O}_7$ solution.

Rinse the burette with (1/10) $\text{K}_2\text{Cr}_2\text{O}_7$ solution and then fill it. 25 ml of the Mohr's salt solution is pipetted out in a 250 ml conical flask, diluted to 150 ml with distilled water, 5 ml of conc. H_2SO_4 , 5 ml of syrupy phosphoric acid are added followed by 4-5 drops of Ba DS indicator. The mixture is titrated with standard $\text{K}_2\text{Cr}_2\text{O}_7$ solution with constant shaking until the colour of the solution just changes from green to reddish-violet.

Table-2 : Titration of Fe^{2+} solution vs standard $\text{K}_2\text{Cr}_2\text{O}_7$ solution

SL. No.	Volume of Fe^{2+} solution taken (ml)	Volume of $\text{K}_2\text{Cr}_2\text{O}_7$ solution required (ml)	Mean volume of $\text{K}_2\text{Cr}_2\text{O}_7$ (ml)
			V ml

Calculation

$$25 \text{ ml of } Fe^{2+} \text{ solution} \equiv V \text{ ml of } f\left(\frac{H}{10}\right) K_2Cr_2O_7 \text{ solution}$$

$$\equiv V \times f \text{ ml of } \left(\frac{H}{10}\right) K_2Cr_2O_7$$

1000 ml of $(\frac{N}{10}) K_2Cr_2O_7$ solution \equiv 1000 ml of $(\frac{N}{10}) Fe^{2+}$ solution
 \equiv 5.585 gm of Fe^{2+}

$\therefore 1 \text{ ml of } (N/10) K_2Cr_2O_7 \equiv 0.005585 \text{ gm of } Fe^{2+}$

$$V \times f \text{ ml of } \left(\frac{N}{T_0}\right) K_2 Cr_2 O_7 \equiv V \times f \times 0.005585 \text{ gm of Fe}^{2+}$$

Now,

25 ml of Fe^{2+} solution contains $V \times f \times 0.005585$ gm Fe^{2+}